

week1

June 15, 2023

Notebook file for Week 1

Import libraries. The nmrbase folder needs to be located within the folder containing this Jupyter notebook

```
[ ]: import numpy as np
import matplotlib.pyplot as plt
import nmrbase.expbases as expbase
import nmrbase.expdata as expdata
```

Tuning response experiment

```
[ ]: filename = r"../DIRECTORY/FILENAME" #Defines the path to the data file. Here
↳the path is relative to the current folder.
```

```
[ ]: a = expbase.expbases() # create an experiment base object (a)
a.load(filename) # load data into a
f1 = plt.figure()
ax1 = f1.subplots() # create axes for a figure

a.plottm(ax1,1) # plot the time domain data of the
↳first scan from the data set
```

```
[ ]: ## ADDITIONAL TASKS:
## include statements such as ax1.set_ylabel() to label axes
## use statements to change appearance such as font size, etc.
## Adjust horizontal axis limits to display only the ring-down signal with
↳statements such as ax1.set_xlim([0,1]). You will need to select appropriate
↳values.

f1 # show the figure again with modifications
```

```
[ ]: f2=plt.figure()
ax2=f2.subplots()

a.pproc['ftmin']=0 # time interval for Fourier transform
↳(0 = all data)
a.pproc['ftmax']=0
```

```

a.pproc['ffmin']=0                                # frequency interval for spectrum
↳display (0 = all data)
a.pproc['ffmax']=0
a.proc()                                           # calculate Fourier transform of the data
↳in a

a.plotfrq(ax2,1)                                  # plot the frequency domain data of the
↳first scan from the data set

## TASKS:
## change parameter "ftmin" and "ftmax" to select the time interval
↳corresponding to the ringdown for the Fourier transform
## change parameter "ffmin" and "ffmax" to select the frequency range of the
↳tuning response signal
## include statements to label axes
## use statements to change appearance such as font size, etc.
## Adjust horizontal axis limits to display the spectrum of the ring-down signal

```

```

[ ]: i=np.argmax(a.frq[0].dta)                    # find data point index of largest peak in
↳spectrum
f=a.frq[0].ind_to_x(i)                            # find frequency corresponding to largest peak in
↳spectrum
print("Tuning frequency: {:.2f} Hz".format(f))     # this statement prints the
↳identified frequency in a nice format

```

```

[ ]: # plot figure again with line indicating maximum peak position. Notice this
↳works only if previous parameters have been correctly selected.
ax2.axvline(x = f, color = 'r', label = 'Larmor Frequency')

f2 # show the figure again with the line included

```

Pulse length calibration. NOTE: To analyze and report multiple datasets, simply copy the necessary blocks of code and change the name of the variables.

```

[ ]: filename2 = r"../DIRECTORY/FILENAME"         # path of data file

```

```

[ ]: a2 = expbase.expbase()
a2.load(filename2)                                # load data

f3=plt.figure()
ax3=f3.subplots()
a2.plottm(ax3,5)

```

```

[ ]: ## TASKS:
## fine tune with ax3.set_xlim and ax3.set_ylim parameters to zoom in on the
↳echo

```

```
f3    # display adapted figure
```

```
[ ]: a2.pproc['digfmin']=1500    # set appropriate digital filter parameters
a2.pproc['digfmax']=3500
a2.digfilt()    # perform digital filter operation
```

```
[ ]: f4=plt.figure()
ax4=f4.subplots()
a2.plottm(ax4,5)    # plot the 5th scan (the 1st scan doesn't have signal
    ↳because the pulse is too short)
```

```
[ ]: ## TASKS:
## fine tune with ax4.set_xlim and ax4.set_ylim parameters to zoom in on the
    ↳echo
```

```
f4    # display adapted figure
```

```
[ ]: f5=plt.figure()
ax5 = f5.subplots()

a2.pproc['ftmin']=0    # time interval for Fourier transform
    ↳(0 = all data)
a2.pproc['ftmax']=0
a2.pproc['ffmin']=0    # frequency interval for spectrum
    ↳display (0 = all data)
a2.pproc['ffmax']=0
a2.proc()    # calculate Fourier transform of the data
    ↳in a
## this will use the digitally filtered data from before. Instead, the original
    ↳data can be processed by loading it again.

a2.plotfrq(ax5,5)    # plot the frequency domain data of the
    ↳first scan from the data set
```

```
## TASKS:
## change parameter "ftmin" and "ftmax" to select the time interval
    ↳corresponding to the ringdown for the Fourier transform
## change parameter "ffmin" and "ffmax" to select the frequency range of the
    ↳tuning response signal
## include statements to label axes
## use statements to change appearance such as font size, etc.
## Adjust horizontal axis limits to display the spectrum of the ring-down signal
```

```
[ ]: # This section is to integrate the peaks in the NMR spectra to determine pulse
    ↳length
f5=plt.figure()
ax5=f5.subplots()
```

```

a2.pproc['intmin']=0      # set correct frequency range for integration
a2.pproc['intmax']=0

a2.integrate()            # perform integration

#find the starting pulselength and increment, then set the correct x-axis
dx=a2.pinc["inc"][0]
x0=a2.p["p1"]
print('x0 =',x0,'s , dx =',dx,'s')      # x0 is the starting pulse length,
    ↳and dx is the increment
a2.idt.x0=x0
a2.idt.dx=dx

a2.idt.plot(ax5,disp=[0])      # disp=0 plots only the NMR signal,
    ↳in trace 0
ax5.axvline(x = 0.001, color = 'r', label = '90° pulse length') # display
    ↳vertical line at 90 degree pulse length

## TASKS:
## change the intmin and intmax parameters to select the correct frequency,
    ↳range for integration
## use set_xlabel and set_ylabel to set the labels to get a publication quality,
    ↳figure
## adjust appearance of figure as needed
## change and report the "x" parameter to indicate the pulse length for 90°,
    ↳pulse with a vertical line

```