

week2b

August 22, 2023

Notebook file for Week 2, T2 Relaxation experiment

Import libraries. The nmrbase folder needs to be located within the folder containing this Jupyter notebook

```
[ ]: import numpy as np
import scipy
import matplotlib.pyplot as plt
import nmrbase.expbases as expbases
import nmrbase.expdata as expdata
```

T2 Relaxation NOTE: To analyze and report multiple datasets, simply copy the necessary blocks of code.

```
[ ]: filename = r"../DIRECTORY/FILENAME"  #Defines the path to the data file. Here
      ↳the path is relative to the current folder.
```

```
[ ]: a = expbases.expbases()
a.load(filename)                                # load data

f1=plt.figure()
ax1=f1.subplots()
a.split()                                       # splitting raw T2 data for later
      ↳processes
a.plottm(ax1,0)
```

```
[ ]: ## TASKS:
## fine tune with set_xlim and set_ylim parameters to zoom in on the echo
ax1.set_ylim([-0.15,0.15])
ax1.set_xlim([0.15,0.25])

f1
```

```
[ ]: a.pproc['digfmin']=1500    # set appropriate digital filter parameters
a.pproc['digfmax']=3500
a.digfilt()                        # perform digital filter operation
```

```
[ ]: ## TASKS:
## fine tune with set_xlim and set_ylim parameters to zoom in on the echo
```

```
f2=plt.figure()
ax2=f2.subplots()

a.plotm(ax2,0)    # plot the entire scan
ax2.set_ylim([-0.15,0.15])
ax2.set_xlim([0.15,0.25])    #zoom in on the 1st echo, which is the strongest one
```

```
[ ]: f3=plt.figure()
ax3 = f3.subplots()

a.pproc['ftmin']=0          # time interval for Fourier transform (0 =
    ↳all data)
a.pproc['ftmax']=0
a.pproc['ffmin']=0          # frequency interval for spectrum display
    ↳(0 = all data)
a.pproc['ffmax']=0
a.proc()                    # calculate Fourier transform of the data
    ↳in a
## this will use the digitally filtered data from before. Instead, the original
    ↳data can be processed by loading it again.

a.plotfrq(ax3,0)            # plot the frequency domain data of the
    ↳first scan from the data set

## TASKS:
## change parameter "ftmin" and "ftmax" to select the time interval
    ↳corresponding to the echoes for the Fourier transform
## change parameter "ffmin" and "ffmax" to select the frequency range of the
    ↳NMR signal
## include statements to label axes
## use statements to change appearance such as font size, etc.
```

```
[ ]: # Time-resolved integration of peaks on NMR spectra
f4=plt.figure()
ax4=f4.subplots()

a.pproc['intmin']=0         # set correct frequency range for integration
a.pproc['intmax']=0

a.integrate()               # perform integration

#find the waiting time tau, then set the correct x-axis
dx=a.pproc["soff"]
x0=a.pproc["slen"]
print('x0 =',x0,'s , dx =',dx,'s')    # x0 is the time of 1st echoes, and
    ↳dx is the interval between echoes
```

```

a.idt.x0=x0
a.idt.dx=dx

a.idt.plot(ax4,disp=[0])           # disp=0 plots only the NMR signal
    ↳ in trace 0

ax4.set_xlabel("seconds [s]")
ax4.set_ylabel("Y LABEL")

## TASKS:
## change the intmin and intmax parameters to select the correct frequency
    ↳ range for integration
## use set_xlabel and set_ylabel to set the labels to get a publication quality
    ↳ figure

```

Fitting of the integral to time axis. NOTE: For some results to fit, optimize the initial conditions

```

[ ]: # Fitting of the integral

ax5 = plt.figure().subplots()
a.idt.plot(ax5,disp=[0])           # to be overlayed by fitted curve,
    ↳ disp=0 plots only the NMR signal in trace 0

def fun(t,a,b,c):
    return FORMULA                 # INPUT the formula used for
    ↳ fitting. Use "np.exp()" for exponential, and t for time-axis.

y = a.idt.dta[0]
x = np.linspace(a.idt.x0,a.idt.x0+a.idt.dx*round(len(y)-1),len(y))

a.p1,a.p2=scipy.optimize.curve_fit(fun,x,y,p0=[y[0],0.1,0],maxfev=5000)
    ↳ # OPTIMIZE the initial conditions

print('T\u2082 = ',round(1/a.p1[1],SF1),'±',round(np.linalg.eig(a.p2)[0][1]**0.
    ↳ 5,SF2),'s')                  # REPLACE "SF1" and "SF2" with positive integers to
    ↳ report the fitting results with correct significant figures

x=np.linspace(x[0],x[-1],1000)     # use 1000 points to
    ↳ generate a smooth curve
a.fit_points = a.p1[0]*np.exp(-a.p1[1]*x)+a.p1[2]
pl,=ax5.plot(x,a.fit_points,'r-')

ax5.set_xlabel("seconds [s]")       # SET the labels to get a
    ↳ publication quality figure
ax5.set_ylabel("Y LABEL")
pl.figure.set_tight_layout('pad')
pl.figure.canvas.draw()

```

```
## TASKS:  
## change the initial fitting parameters to obtain proper fitting  
## report the fitted relaxation constants with proper significant figures  
## use set_xlabel and set_ylabel to set the labels to get a publication quality  
↪ figure  
## adjust appearance of figure as needed
```